

1 A novel analytical pyrolysis device applicable for measurements of less volatile pyrolyzates

2

3 Hirohisa Kano^a, **Takahiro Okamoto^a**, Shinya Kitagawa^a, Yoshinori Iiguni^a, Hajime Ohtani^{a,*}, Hiroshi Ito^b,

4 Koichiro Iwai^{a,b}, Minoru Kuno^c

5

6 ^a *Department of Materials Science and Engineering, Graduate School of Engineering, Nagoya Institute of*

7 *Technology, Nagoya 466-8555, Japan*

8 ^b *Toyota Central R&D Labs., Inc., Nagakute, Aichi, 480-1192, Japan*

9 ^c *GL Science Inc., 237-2, Sayamagahara, Iruma, Saitama, 358-0032, Japan*

10

11 ABSTRACT

12 A novel analytical pyrolysis device suitable for the measurements of the less volatile pyrolyzates was

13 developed. This device is comprised of a pyrolysis unit based on the vertical microfurnace pyrolyzer and a

14 movable pyrolysis tube made of stainless steel. The upper side of the pyrolysis tube is first set at the heated

15 center of the pyrolysis unit. A polymer sample in a sample cup is introduced into the heated center and

16 pyrolyzed under the flow of nitrogen carrier gas. The pyrolyzates are then transferred to the bottom side of

17 the pyrolysis tube which is cooled down by air blow and the less volatile pyrolyzates are condensed. After

18 that the pyrolysis tube is moved down out of the pyrolysis unit and blocked from the hot zone by moving

19 the three port valve. Finally the trapped pyrolyzates **are** flushed by an appropriate solvent into a collection

20 vial. The recovered pyrolyzates in the solution can be then subjected to the measurements suitable for less

21 volatile products such as matrix-assisted laser desorption/ionization mass spectrometry. Moreover, this

22 **device** is also on-line connected to the injection port of a high performance liquid chromatographic system,

23 which is also appropriate for the measurements of oligomeric products. The basic performance of the

24 developed **device** was tested by using some polymer samples forming less volatile pyrolyzates such as

25 **poly(methyl acrylate), Nylon 11,** poly(butylene terephthalate) and poly(2,6-dimethyl-1,4-phenylene ether).

26 **In addition, the usefulness of this method was demonstrated through the characterization of the crosslinking**

27 sequences in a photocured resin combined with the thermally-assisted hydrolysis and methylation reaction.

28

29 *Keywords:* Less volatile pyrolyzate, HPLC, MALDI-MS, Pyrolysis device, Analytical pyrolysis, Polymer
30 characterization

31

32 1. Introduction

33 Recently, modern pPyrolysis-gas chromatography (Py-GC) is widely utilized ~~for~~ in the field of
34 polymer characterization not only for simple “fingerprint” identification and compositional analysis, but
35 also for the detailed structural analysis such as copolymer sequences, region- and stereoregularities, end
36 groups, and branching and cross-linking structures of various polymeric materials [1-3]. For such advanced
37 objectives, it is often demanded to inspect a wide range of pyrolyzates including large molar mass
38 (oligomeric) products and/or highly polar compounds. This technique, however, is inherently inapplicable
39 for the measurements of such less volatile pyrolyzates ~~such as large molar mass (oligomeric) products and~~
40 ~~highly polar compounds~~ due to the use of GC separation, even if they possess useful information for the
41 characterization of the sample materials. Therefore, an analytical pyrolysis technique suitable for such
42 intractable pyrolyzates has been required.

43 A possible approach to overcome this issue ~~should be~~ is the combination of pyrolysis and high
44 performance liquid chromatographic (HPLC) separation ~~applicable even for less volatile products~~.
45 Stepwise (off-line) analyses of collected pyrolyzates by solvent-trap or extraction have been reported by
46 using HPLC (or supercritical fluid chromatography (SFC)) for polymeric samples such as polystyrene [4,
47 5] and latex containing polar vinyl acid components [6]. In 1996 van der Hage and Boon [7] proposed an
48 on-line Curie-point pyrolysis-HPLC-mass spectrometry (MS) system. In this ~~system~~ method, by controlled
49 solvent flow with appropriate valve switching, condensed pyrolyzates in the pyrolysis chamber were
50 washed out into a precolumn for trapping. The collected products were then flushed into the separation
51 column of HPLC by switching the solvent and its flowing path with the valve control. The feasibility and
52 the potential of this ~~system~~ technique were demonstrated through the analysis of high molar mass pyrolysis

1 53 products of lignin polymers. Recently, a generic pyrolysis interface called “PyroVial” was presented [8].
2
3 54 The sample in the closed vial was electrically heated up to 1000 °C on the dedicated socket. Combined
4
5 55 with a sophisticated syringe manipulation, the vial containing the pyrolyzates can be hyphenated to any
6
7 56 chromatographic methods including GC, HPLC and SFC.

8
9
10 57 Alternative measure for this subject would be the direct combination between pyrolysis and MS
11
12 58 without separation. In the conventional pyrolysis-MS technique, however, vaporized pyrolyzates have been
13
14 59 predominantly measured using electron ionization (EI), chemical ionization (CI) and field ionization (FI)
15
16 60 sources [9]. Less volatile pyrolyzates, therefore, cannot be analyzed and considerable fragmentation
17
18 61 sometimes makes the observed mass spectra too complicate to interpret especially for EI source.
19
20 62 Meanwhile, matrix assisted laser desorption/ionization (MALDI)-MS and electrospray ionization (ESI)-MS
21
22 63 are widely utilized in recent years for the characterization of ~~oligomeric and polymeric~~ large molar mass
23
24 64 compounds especially for oligomeric components in recent years [10, 11]. Thus the combination between
25
26 65 pyrolysis and MALDI- and/or ESI-MS ~~must be~~ is also a promising approach to obtain the information of
27
28 66 less volatile pyrolyzates.

29
30
31
32
33
34 67 In this work, a novel pyrolysis device suitable for the analysis of the less volatile pyrolyzates was
35
36 68 developed on the basis of a microfurnace pyrolyzer which is one of the representative pyrolyzers used for
37
38 69 analytical purpose [3]. Pyrolysis reaction is achieved basically under the same conditions for conventional
39
40 70 Py-GC measurements. The device can trap the less volatile pyrolyzates in a movable pyrolysis tube and
41
42 71 flush the trapped components with an appropriate solvent. The recovered pyrolyzates in the solvent are
43
44 72 readily analyzed by MALDI- and ESI-MS. Furthermore, this device can be also on-line connected to the
45
46 73 injection port of HPLC. The basic performance and applicability of the device was tested by using some
47
48 74 polymer samples forming less volatile pyrolyzates.

49
50
51
52
53 75

54 76 **2. Experimental**

55 77 *2.1. Construction and operation of pyrolysis device*

56
57 78 Fig. 1 shows the schematic diagram and the operation of the developed pyrolysis device. The

~~operation of the device was shown in Fig. 2.~~ This device is comprised of a pyrolysis unit based on a vertical microfurnace pyrolyzer [3], of which temperature is controlled by a PID controller (EC5530, Ohkura) and a movable pyrolysis tube (5) made of stainless steel (190 mm long x 4 mm o.d. x 2 mm i.d.) set inside of an outer tube (1/4 inch o.d. x 4.35 mm i.d.). The upper side of the pyrolysis tube is first positioned at the heated center of the pyrolysis unit (Fig. 2(a)). A polymer sample weighed about 1 mg in a surface deactivated stainless steel sample cup (6) (Eco cup, Frontier Lab) is introduced into the heated center at 500 or 550 °C and pyrolyzed under the flow of nitrogen carrier gas (70 mL min⁻¹). ~~As for thermally-assisted hydrolysis and methylation (THM) of a photocured resin sample, 3 mg of the resin mixed with 5 μL of a 25 wt% of tetramethylammonium hydroxide (TMAH) (Aldrich) in the sample cup was introduced into the pyrolysis unit set at 400 °C.~~ The pyrolyzates are then transferred to the bottom side of the pyrolysis tube which is cooled down by cooling fans (4) and the less volatile pyrolyzates are condensed. To increase the trapping efficiency, the connecting union at the bottom of the pyrolysis tube is packed with quartz wool (10) (3 mg). During this step nitrogen seal gas (10 mL min⁻¹) (3) flows to prevent the leakage of pyrolysis products outside of the pyrolysis (inner) tube. After that the pyrolysis tube is moved down out of the pyrolysis unit sealed with a tough O-ring (9) (Kalrez S004, 3.5 mm x 1.5 mm, DuPont) and blocked from the hot zone by switching the three port valve (8) (Fig. 2(b)).

Fig. 1. Schematic diagram and operation of pyrolysis unit with movable trapping tube for pyrolyzates.

Fig. 2. Operation of the pyrolysis unit.

Finally the trapped pyrolyzates are flushed by an appropriate solvent (7) into a collection vial (5 mL volume) (11). The flushing solvent (1 mL THF; HPLC grade; Tokyo Chemical Industry) is pushed out to the pyrolysis tube with the seal gas by rotating another 3-port switching valve (8) for the seal gas. The recovered pyrolyzates in the solution ~~can be~~ are readily subjected to the measurements for less volatile products such as MALDI- and ESI-MS. Moreover, this ~~system device~~ can be also on-line connected to the injection port of HPLC. After the PTFE 3-port valve (15) is switched and the position of the collection vial

1 105 is adjusted, the solution of the recovered pyrolyzates is transferred to the sample loop (2 μ L volume) of the
2
3 106 injector (Rheodyne 7725i) with the pressure of the seal gas. ~~(A photograph of the system and the detailed~~
4
5 107 ~~operation are given as supplementary data.)~~
6

7 108

10 109 2.2. Conditions for on-line HPLC

11
12 110 On-line HPLC measurements of the collected pyrolyzates were carried out with an LC-VP HPLC
13
14 111 system (Shimadzu) consisting of dual solvent delivery pumps (LC-10ADvp; Shimadzu), a degasser
15
16 112 (DGU-12A, Shimadzu), a column oven (CO-10ADvp, Shimadzu), a photodiode array UV-visible detector
17
18 113 (SPD-M10A, Shimadzu), and a system controller (SCL-10A, Shimadzu). The separation column was
19
20 114 Chemcobond 3.5-ODS-H (Chemco, 150 mm x 1.0 mm, 3.5 μ m particles) connected with a guard column
21
22 115 of the same packings (Chemco, 35 mm x 1.0 mm). A mixed solvent of water and THF was adopted as the
23
24 116 mobile phase for a gradient elution. The gradient conditions were programmed from 30 (for PBT) or 20%
25
26 117 (for PPE) THF (10 min hold) to 100% THF (for 60 min), occasionally added with 0.1% (v/v) trifluoroacetic
27
28 118 acid mainly for polar products from PBT. The column oven was set at 40 $^{\circ}$ C. The UV detection was carried
29
30 119 out either at 240 (for PBT) or 280 (for PPE) nm.
31
32
33
34
35

36 120 2.3. Conditions for MALDI-MS

37
38 121 ~~MALDI-MS measurements of the collected pyrolyzates were performed to confirm the high molar~~
39
40 122 ~~mass oligomeric and polar products that cannot be detected by Py-GC.~~ In this work, either high resolution
41
42 123 spiral time-of-flight (TOF)-MS (JMS-S3000, JEOL) or conventional TOF-MS (AXIMA-CFRplus,
43
44 124 Shimadzu Kratos) was used. The spiral TOF-MS was equipped with Nd:YLF laser (349 nm) and performed
45
46 125 with 0.5 ns pulse width and 250 Hz frequency, while the conventional TOF-MS was with nitrogen laser
47
48 126 (337 nm), operated in 3 ns pulse width and 10 Hz frequency. Ions generated by the laser desorption were
49
50 127 introduced into the flight tube with an accelerating voltage of 20 kV and separated in a positive spiral mode
51
52 128 (JMS-S3000) and a linear mode (AXIMA-CFR plus). The mass spectra were acquired by averaging 30-100
53
54 129 and 100 individual laser shots, with the delay time set at 150 and 80 ns, respectively for the spiral and the
55
56
57
58
59
60 130

1 131 conventional TOF. Mass calibration was accomplished either by using poly(ethylene glycol) (average
2
3 132 molecular weight 400 - 2000) as external standard or a hindered amine light stabilizer (Chimassorb 119FL,
4
5 133 Ciba Specialty Chemicals, m/z 2285.61 for the monoisotopic ion) as the internal standard. As for the matrix
6
7 134 reagent, 2-(4-hydroxyphenylazobenzoic acid) (HABA; Sigma-Aldrich), 4-hydroxybenzylidene
8
9
10 135 malononitrile (HBMN; Alfa Aesar) and
11
12 136 trans-2-[3-(4-tert-butylphenyl)-2-methyl-2-propenylidene]malononitrile (DCTB; Tokyo Chemical Industry)
13
14
15 137 were used. Each matrix solution (10 mg mL^{-1}) was prepared by dissolving 1 mg of the matrix reagent in 0.1
16
17 138 mL of THF (HPLC grade; Tokyo Chemical Industry), then 10 μL each of the matrix and the sample
18
19 139 solutions were mixed, and an aliquot (1 μL) of the mixed solution dropped on the sample plate was
20
21
22 140 subjected to MALDI-MS measurement. In the case of internal calibration, 1 μL of the calibrant solution (1
23
24 141 mg mL^{-1} in THF) was also added to the sample/matrix mixed solution. The ions observed in the MALDI
25
26 142 mass spectra were predominantly sodium-cationized molecules of the products under the conditions in this
27
28
29 143 work.

30 31 144 32 33 34 145 *2.4. Conditions for HPLC-ESI-MS*

35
36 146 In order to identify the collected pyrolyzates of poly(phenylene ether) (PPE) subjected to HPLC
37
38 147 measurements, a HPLC-MS system consisting of UPLC system (Acuity, Waters) and TOF-MS (SYNAPT
39
40
41 148 G2, Waters) was used. MS was operated in the positive ESI mode under the conditions of source
42
43 149 temperature 120 $^{\circ}\text{C}$, desolvation N_2 gas temperature 350 $^{\circ}\text{C}$, desolvation gas flow 800 L h^{-1} , capillary
44
45 150 voltage of 3.0 kV, and cone voltage 40 V. The separation column and the conditions for measurements were
46
47
48 151 basically the same as those for the on-line HPLC-UV detection (section 2.2).

49 50 152 51 52 53 153 *2.5. Conditions for Py-GC-MS*

54
55 154 Py-GC-MS measurements of the polymer samples were also performed for comparison. A vertical
56
57 155 microfurnace pyrolyzer (PY-2020iD, Frontier Lab) was attached to the injection port of a GC-MS system
58
59
60 156 (GCMS-QP2010 Plus, Shimadzu). About 0.1 mg of a polymer sample weighed in a deactivated stainless

1 157 steel sample cup was pyrolyzed at 500 or 550 °C under the flow of helium carrier gas. The flow rate of 50
2
3 158 mL min⁻¹ of carrier gas at pyrolyzer was reduced to 1.0 mL min⁻¹ at a metal capillary separation column
4
5 159 (Ultra-Alloy 5 MS, 30 m x 0.25 mm i.d., coated with 0.25 μm of 5% diphenyl-95% dimethyl-polysiloxane,
6
7 160 Frontier Lab) by means of a splitter. The temperature both of PY/GC interface and the GC injection port
8
9
10 161 was set at 320 °C. The column oven temperature was first held at 40 °C for 2 min, then programmed to
11
12 162 320 °C at a rate of 10 °C min⁻¹ and finally held at 320 °C for 30 min. The separated components were
13
14
15 163 introduced into the ionization source of a quadrupole MS through a GC/MS interface held at 320 °C. The
16
17 164 ionization was achieved with electron ionization at 70 eV. The mass scan range was *m/z* 29 - 600.

19 165 20 21 22 166 *2.6. Polymer samples*

23
24 167 Poly(methyl acrylate) (PMA) sample was radically polymerized in a standard procedure. Industrially
25
26 168 available pelletized samples were used for Nylon 11 (polyundecanoamide) and poly(butylene terephthalate)
27
28
29 169 (PBT). Poly(2,6-dimethyl-1,4-phenylene ether) (PPE) sample obtained using Mn-based catalyst was the
30
31 170 same as that used in previous work [12]. The results for PMA and Nylon 11 are discussed as supplementary
32
33
34 171 information. ~~The chemical structure of the polymer samples are shown in Fig. 3.~~

35
36 172 ~~The photocured acrylic resin sample was prepared with neopentyl glycol diacrylate through~~
37
38 173 ~~photopolymerization in the presence of 7 wt% 2,2'-dimethoxy-2-phenylacetophenone as a photoinitiator~~
39
40
41 174 ~~using the UV-cured equipment consisting of a conveyor and a medium-pressure mercury lamp.~~

42
43 175
44
45 176 **Fig. 3. Polymer samples.**

46 47 177 48 49 178 **3. Results and discussion**

50 51 52 179 *3.1. Poly(methyl acrylate) (PMA)*

53
54 180 ~~—The basic performance of the pyrolysis device was tested for PMA sample. PMA generally~~
55
56 181 ~~decomposes into a series of MA oligomers through the random chain scission accompanied by radical~~
57
58
59 182 ~~transfer reactions at elevated temperature [13, 14]. Fig. 4 shows a typical pyrogram of PMA observed by~~

183 conventional Py-GC-MS at 550 °C. In this case MA monomer and oligomers up to heptamer ($n = 7$) can be
184 observed [3]. Recently, by using a special Py high temperature GC equipped with a septum free injector
185 [15] and a thermally stable short capillary column coated with carborane siloxane stationary phase, it was
186 reported that the larger oligomers up to 16-mer were detected in the pyrogram of PMA under the conditions
187 for the GC system around 400 °C [16].

Fig. 4. Pyrogram of PMA observed by Py-GC-MS at 550 °C.

188
189
190
191 — PMA sample was then subjected to the developed pyrolysis device. An aliquot of the collected
192 pyrolyzates (at 550 °C) dissolved in THF was measured by MALDI spiral-TOF-MS using HABA as a
193 matrix reagent. The observed mass spectrum is shown in Fig. 5. A series of MA oligomers are observed at
194 least up to 26-mer, and each n -mer contains three types of products as shown in Fig. 5(b). These results
195 demonstrate that high molar mass pyrolyzates unobservable by Py-GC are successfully recovered by the
196 pyrolysis device.

Fig. 5. MALDI mass spectrum of pyrolysis products at 550 °C of PMA collected in THF

3.2. Nylon 11

197
198
199
200
201 — Aliphatic polyamides (Nylon) are known to be pyrolyzed mainly into various sizes of cyclic amides
202 (lactams) and linear products having nitrile, alkyl and/or olefinic terminals formed through *cis*-elimination
203 reactions along with random chain scissions. As for Nylon 11, therefore, the pyrolysis products listed in
204 Table 1 should be yielded. Fig. 6 shows a typical pyrogram of Nylon 11 obtained by Py-GC-MS at 550 °C.
205 Although various pyrolyzates listed in Table 1 are observed [3, 17], the assigned products are mostly those
206 having only one or no amide group such as L_1 , L_2 and D_{0-10} . Even the largest pyrolyzate, A_{1-20} , of which
207 molecular weight is *ca.* 508, contains merely two amide groups.—

209 **Table 1** Representative pyrolyzates of Nylon 11.

210 **Fig. 6.** Pyrogram of Nylon11 observed by Py-GC-MS at 550 °C.

211
212 ~~— Nylon 11 was then pyrolyzed at 550 °C in the pyrolysis device and the pyrolyzates were collected with~~
213 ~~THF. The recovered solution was sampled to be measured directly by MALDI MS. Fig. 7 shows the~~
214 ~~observed mass spectra for the pyrolyzates of Nylon 11 obtained by MALDI spiral TOF MS with HABA~~
215 ~~matrix. In the oligomer range larger than m/z 500, which cannot be detected by Py-GC MS, a number of~~
216 ~~products are detected. Among these, trimeric and tetrameric lactams (L_3 and L_4) are distinctly observed. As~~
217 ~~shown in the expanded spectrum in Fig. 7(b), various linear products containing three amide groups were~~
218 ~~identified in the m/z range between 620 and 645. (The observed peak splitting should be caused by the~~
219 ~~overlapping of possible isotopic ions for individual pyrolyzates.) Moreover, although with small intensities,~~
220 ~~ions for the pyrolyzates were clearly observed even in the region around 1400, which corresponding to the~~
221 ~~octameric products with seven amide groups. These observations indicated that not only large molar mass~~
222 ~~but also polar pyrolyzates originated from polar polymers such as Nylons should be able to be analyzed by~~
223 ~~using the developed pyrolysis device.~~

224
225 **Fig. 7.** MALDI mass spectrum of pyrolysis products at 550 °C of Nylon11 collected in THF.

226 *3.1 Basic performance of the developed pyrolysis device*

227 *3.1.1. Poly(butylene terephthalate) (PBT)*

228 The basic performance of the developed pyrolysis device including on-line HPLC measurements of
229 the pyrolyzates obtained by this device ~~was then performed~~ evaluated using PBT, which is widely utilized
230 as one of typical engineering plastics and contains aromatic moieties suitable for UV detection. In general,
231 polyesters like PBT thermally decompose into oligomeric products primarily through *cis*-elimination at the
232 ester linkages followed by partial decarboxylation. It has been reported for PBT, therefore, several types of
233 products listed in Table 1 should be formed by pyrolysis reactions ~~initiated chiefly with the *cis* elimination~~
234 ~~shown in Fig. 8 [12-14].~~

235

236

Table 1 Possible pyrolyzates of PBT.

237

Fig. 8. Pyrolysis pathway of PBT.

238

239 Fig. 2 shows a pyrogram of PBT obtained by Py-GC-MS at 500 °C. The main pyrolyzates observed

240 are 1,3-butadiene and A_n, B_n and C_n with degree of polymerization (n) of 0 or 1 [3, 12]. The yields of E_n

241 and F_n products having phenyl terminal(s) were fairly small, probably because the contribution of

242 decarboxylation to form phenyl terminals would be minor under the conditions of relatively low pyrolysis

243 temperature at around 500 °C. Meanwhile highly polar decarboxylic acid products (D_n) only shows a broad

244 hump around 20 min for terephthalic acid (D₀). At least the larger pyrolyzates than C₁, or the products with

245 molar mass > 500, were not able to be observed by conventional Py-GC-MS at all.

246

247 **Fig. 2.** Pyrogram of PBT observed by Py-GC-MS at 500 °C.

248

249 Next PBT was pyrolyzed at 500 °C in the pyrolysis device, and the pyrolyzates were collected in THF.

250 The obtained THF solution was first measured directly by conventional MALDI-MS with HBMN as a

251 matrix reagent. Fig. 3 shows the MALDI mass spectrum of the recovered pyrolyzates of PBT. The

252 pyrolyzates observed in the pyrogram by Py-GC-MS (Fig. 2) such as B₀, A₁-C₁ and E₁ were also observed,

253 although smaller products (A₀, E₀, etc) might evaporate during MALDI process and not be detected in the

254 mass spectrum. In addition, the pyrolyzates that could be hardly detected by Py-GC-MS, *e. g.* D₁, F₁ and

255 the products with greater degrees of polymerization (n = 2, 3), were also clearly confirmed in the MALDI

256 mass spectrum (Fig. 3). These observations prove that such larger pyrolyzates from PBT can be sufficiently

257 recovered by the pyrolysis device.

258

259 **Fig. 3.** MALDI mass spectrum of pyrolysis products at 500 °C of PBT collected with THF.

260

261 The solution of the collected pyrolyzates was then subjected to on-line HPLC measurements. Fig. 4

262 shows the UV-detected chromatogram of the recovered pyrolyzates of PBT. In this chromatogram, several
263 peaks were clearly observed. In order to identify the peak components, the eluents around the main peaks
264 (26 - 29, 34 - 37, 42 - 44, and 46 - 48 min) were fractionated (5 times for accumulation) and subjected to
265 MALDI-MS measurements. As the results, the pyrolysis products designated in Fig. 4 were at least
266 assigned. Among these, the highly polar and/or large molar mass components, which cannot be detected by
267 Py-GC-MS, were included (A₂, B₂, C₂, D₂, E₂, etc). This result demonstrates that on-line HPLC
268 measurements of less volatile pyrolyzates are able to be performed by using the analytical pyrolysis system
269 developed.

Fig. 4. On-line Py-HPLC-UV chromatogram of collected pyrolysis products of PBT at 500 °C.

3.1.2. Poly(2,6-dimethyl-1,4-phenylene ether) (PPE)

Finally The pyrolysis device was **then** applied to the characterization of PPE. It has been reported that
PPE thermally decompose into various oligomers in a unique manner with the major contribution of Fries
rearrangement reaction to form C-C bondings along with the simple cleavages of the main chains as shown
in Fig. 5 [15, 16]. Fig. 6 shows a typical pyrogram of PPE obtained by Py-GC-MS at 550 °C. Although a
wide range of oligomeric products should be formed, observed pyrolyzates by Py-GC-MS were limited to
mainly monomers, dimers and trimers together with minor amounts of tetramers. Table 2 summarizes the
peak assignments of major pyrolyzates in dimeric and trimeric regions. Reflecting the considerable
contribution of Fries arrangement, various types of products containing some isomers with methylene
and/or ether linkage(s) can be formed [3, 15, 16].

Fig. 5. Typical pyrolysis pathway of PPE .

Fig. 6. Pyrogram of PPE observed by Py-GC-MS at 550 °C.

Table 2 Peak assignment of PPE pyrolyzates observed in dimer and trimer regions in the pyrogram.

The PPE sample was also pyrolyzed in the new pyrolysis devices at 550 °C. The recovered solution **is**

289 ~~THF~~ of the pyrolyzates in THF was first analyzed by MALDI-spiral-TOF-MS using DCTB as the matrix
290 reagent. The observed mass spectrum is shown in Fig. 7. In addition to the pyrolyzates up to tetramers ($n =$
291 4) observed by Py-GC-MS, higher molar mass oligomers are clearly detected at least up to heptamers ($n =$
292 7). This fact shows the less volatile pyrolyzates of PPE can be successfully recovered by the pyrolysis
293 device.

294
295 **Fig. 7.** MALDI mass spectrum of pyrolysis products at 550 °C of PPE collected in THF.

296
297 ~~Here it should be noted that a number of clusters were observed for oligomers of each degree of~~
298 ~~polymerization (n). Hence the peaks in the pentamer region ($n = 5$) were examined in detail. A~~
299 ~~representative formation pathway of a pentameric ($n = 5$) product is shown in Fig. 15, provided that the~~
300 ~~Fries rearrangement would completely proceed during pyrolysis pathway. The resultant pyrolyzate shown~~
301 ~~in Fig. 15 should be possibly the largest pentamer, and is abbreviated to 5-0. Depending on the cleavage~~
302 ~~position, the smaller pentamers by one or two CH_2 unit(s) could be also produced, which are abbreviated~~
303 ~~5-1 and 5-2, respectively. These three types of products, 5-0—5-2, would be thus the most probable~~
304 ~~pentameric pyrolyzates. In fact, as for the dimeric and trimeric pyrolyzates of PPE observed by Py-GC-MS,~~
305 ~~the corresponding products were confirmed to be predominant as shown in Fig. 13 and Table 3. In the~~
306 ~~observed MALDI mass spectrum of the PPE pyrolyzates, however, 5-0—5-2 were not the main products in~~
307 ~~pentamer region as shown in the expanded view in Fig. 14. Smaller components corresponding to 5-3—5-6~~
308 ~~appeared rather abundantly. Moreover, the larger m/z ions by 2 for 5-0—5-6, which might be attributed to a~~
309 ~~hydrogenated ($+ \text{H}_2$) or an oxidized ($+ \text{O} - \text{CH}_2$) product of 5-0—5-6, were detected in higher intensities than~~
310 ~~those for 5-0—5-6. These observations suggest that the undesirable degradations and/or secondary reactions~~
311 ~~might fairly take place for PPE pyrolyzates during MALDI process [23].~~

312
313 **Fig. 15.** Pyrolysis pathway of PPE to form pentamer ($n = 5$) with maximum numbers of substitution (5-0).

1 315 The recovered solution of the PPE pyrolyzates was then subjected to on-line HPLC measurements. Fig.
2
3 316 8 shows the chromatogram of the collected PPE pyrolyzates observed in UV detection at 280 nm. The fact
4
5 317 that a number of separated peaks were clearly observed suggest that the analysis of the less volatile
6
7 318 pyrolyzates for PPE can be achieved by on-line HPLC measurement. In the case of the PPE pyrolyzates,
8
9
10 319 however, the correct identification of each eluent should be difficult by the fractionation-MALDI-MS
11
12 320 technique because of possible decomposition and secondary reaction during laser irradiation as mentioned
13
14
15 321 in supplementary information. Therefore, HPLC-ESI-MS measurements were carried out for the
16
17 322 appropriate identification of the PPE pyrolyzates, in which an aliquot of the collected PPE pyrolyzate
18
19 323 solution was sampled and introduced separately (off-line) into the HPLC-ESI-MS system.
20
21

22 324
23
24 325 **Fig. 8.** On-line HPLC chromatograms of collected pyrolysis products of PPE at 550 °C by UV detection.
25
26 326

27
28 327 ~~Fig. 17 shows the observed total ion chromatogram (TIC) of the collected pyrolyzates of PPE. Since~~
29
30 328 ~~the separation conditions have not been optimized, the observed peak separation is far from sufficient. The~~
31
32
33 329 ~~eluted components, however, can be identified on the basis of the accumulated ESI mass spectrum for the~~
34
35 330 ~~whole range of the chromatogram. Fig. 18 shows the expanded region for the pentamers in the accumulated~~
36
37
38 331 ~~ESI mass spectrum in the range between 10 and 50 min. It is interesting to note that the observed main~~
39
40 332 ~~pentamers (5-0—5-2) are those expected to form by pyrolysis as shown in Fig. 15. This result indicates~~
41
42 333 ~~even the PPE pyrolyzates can be observed in the original form by ESI-MS, unlike MALDI-MS. Moreover,~~
43
44 334 ~~the fact that the ions for pentamers are clearly detected suggests the less volatile pyrolyzates of PPE should~~
45
46 335 ~~be analyzed by HPLC combined with the developed pyrolysis device.—~~
47
48
49

50 336
51 337 **Fig. 17.** Total ion chromatograms of collected pyrolysis products of PPE obtained by HPLC-ESI-MS.
52

53 338 **Fig. 18.** Accumulated ESI mass spectrum of collected pyrolysis products of PPE at 550 °C.
54
55

56 339
57
58 340 ~~Finally,~~ The elution behaviors of individual pyrolyzates were examined in detail based on the
59
60
61
62
63
64
65

1 341 corresponding selected ion chromatograms (SIMs) obtained by HPLC-ESI-MS. Fig. 9 shows the SIMs of
2
3 342 some PPE pyrolyzates for the largest components of individual oligomers (tetramers to hexamers). Not
4
5 343 only the elution times of the products shifted depending on their molecular size, but also the number of the
6
7 344 observed peaks for the individual oligomers increased with increase in m/z . The latter multiple peaks
8
9 345 should be attributed to the possible formation of the isomeric pyrolyzates. As for pentamers (m/z 639), at
10
11 346 least five separated peaks are observed. If the pyrolysis of PPE proceeds with the complete Fries
12
13 347 rearrangement, the formed pentamer ~~(5-0)~~ must be that containing 5 OH groups ~~alone as shown in~~ (Fig. S5).
14
15 348 However, considerable amounts of dimeric and trimeric pyrolyzates containing ether linkage(s) were also
16
17 349 confirmed in the pyrogram of PPE obtained by Py-GC-MS as shown in Fig. 6 and Table 2. This fact
18
19 350 suggests that various isomeric products containing ether linkage(s) should be also formed in the pentamer
20
21 351 region. Even for the pentameric products with the same ~~molecular weight as (5-0)~~ (m/z 639), various
22
23 352 isomeric forms can be considered as shown in Table 3. (Further isomeric products should be produced
24
25 353 depending on the difference in the distributions of OH groups.) Since these pentamers contain various
26
27 354 numbers of OH groups (1 - 5), their polarity should be different from each other which in turn causes the
28
29 355 peak separation in the chromatogram. This result demonstrated that even the isomeric products in less
30
31 356 volatile pyrolyzates can be separated by HPLC measurements which could not be differentiated by MS
32
33 357 analyses alone.
34
35
36
37
38
39
40
41
42

43 359 **Fig. 9.** Extracted ion chromatograms of oligomeric pyrolysis products of PPE at 550 °C.

44 360 **Table 3** Possible isomeric pyrolyzates for pentamer (5-0) of PPE observed in extracted ion chromatogram.

47 361 *3.2 Application to elucidate crosslinking sequences in photocured resin using THM*

48 362 Photocurable acrylic esters are utilized for many industrial materials. Characterization of the network
49 363 structures in the cured resins has been requested to correlate it with various properties of the resins. It has
50 364 been demonstrated that THM-GC often brings us valuable clues on the chemical structure in the photocured
51 365 acrylic resins. In this technique, the ester linkages in the cured resin samples selectively decomposed into
52 366

1 367 methyl derivatives which provide various information such as conversion, copolymer composition and
2
3 368 average molar mass of prepolymers [2]. As for the cross-linking sequences, it was reported that the chain
4
5 369 length distribution of the network junctions composed of up to six acryloyl units in the photocured acrylic
6
7 370 resin sample could be estimated from the peak intensities of the characteristic methyl acrylate (MA)
8
9
10 371 oligomers (up to hexamer) reflecting cross-linking sequence structure [17]. However, the network junctions
11
12 372 containing longer sequence comprised of more than seven acryloyl units were not elucidated because the
13
14
15 373 corresponding higher MA oligomers were not observed in the pyrogram due to their low volatility.

16
17 374 In this case, the developed pyrolysis device could provide the information of the longer cross-linking
18
19 375 sequences in the cured acrylic resins through the THM reaction in the presence of TMAH. Fig. 10 presents
20
21
22 376 possible THM reactions in the photocured resin sample formed from neopentylglycol diacrylate, in which
23
24 377 two types of MA oligomers might be produced. The cured resin sample mixed with TMAH was introduced
25
26
27 378 into the pyrolysis device at 400 °C. The collected products in THF were then analyzed by conventional
28
29 379 MALDI-MS with HABA as a matrix reagent and NaI as a cationization agent. Fig. 11 shows the MALDI
30
31 380 mass spectrum of the THM products of the cured acrylic resin sample. Two series of the MA oligomers
32
33
34 381 were observed at least up to 18-mers. This result demonstrates that the developed pyrolysis device can
35
36 382 provide us the information of the chain length distribution of the network junctions in cured resins even for
37
38
39 383 the longer sequences which cannot be evaluated by THM-GC measurements.

40
41 384
42
43 385 **Fig. 10.** Thermally-assisted hydrolysis and methylation pathway of photo cured acrylic resin sample.

44
45 386 **Fig. 11.** MALDI mass spectrum of THM products of photocured acrylic resin sample collected in THF.

48 387 49 50 388 51 52 53 389 **4. Conclusion**

54
55 390 A novel pyrolysis device suitable for analysis of the less volatile pyrolyzates was developed using a
56
57
58 391 microfurnace type pyrolysis unit combined with a movable pyrolysis tube (made of stainless steel). The
59
60 392 less volatile pyrolyzates of polymer samples were trapped on the inner wall of the pyrolysis tube by air

1 393 cooling, and then solvent-flushed into a collection vial. The recovered solution containing less volatile
2
3 394 pyrolyzates can be directly analyzed by MALDI-MS. The on-line HPLC analyses of the solutions were also
4
5 395 possible. In the present state, **however**, only the HPLC-UV detection system has been on-line connected
6
7 396 with the pyrolysis device. The **more** powerful combination, ~~however~~, could be that between pyrolysis and
8
9
10 397 HPLC-MS. The **development of** on-line connection of ~~a~~-HPLC-ESI-MS ~~system~~ with the pyrolysis device is
11
12 398 currently in progress along with the optimization of the operating conditions.
13
14
15 399

16 17 400 **Acknowledgments**

18
19 401 This work was supported in part by the research foundation of Frontier Laboratory Inc., and **JPSJ**
20
21
22 402 **KAKENHI Grant Number 233350032**.
23

24 403 25 26 404 **References**

- 27
28
29 405 [1] S. C. Moldveanu, Analytical Pyrolysis of Synthetic Organic Polymers, Elsevier, 2005.
30
31 406 [2] T. P. Wampler (Ed.), Applied Pyrolysis Handbook, second ed., CRC Press, 2007.
32
33
34 407 [3] S. Tsuge, H. Ohtani, C. Watanabe, Pyrolysis-GC/MS Data Book of Synthetic Polymers -Pyrograms,
35
36 408 Thermograms and MS of Pyrolyzates-, Elsevier, 2011.
37
38 409 [4] S.-T. Lai, D.C. Locke, J. Chromatogr. A, 255 (1983) 511-527.
39
40
41 410 [5] T. Nonobe, H. Ohtani, T. Usami, T. Mori, H. Fukumori, Y. Hirata, S. Tsuge, J. Anal. Appl. Pyrolysis, 33
42
43 411 (1995) 121-138.
44
45 412 [6] F. C.-Y. Wang, B. Gerhart, C. G. Smith, Anal. Chem., 67 (1995) 3681-3686.
46
47
48 413 [7] E.R.E. van der Hage, J.J. Boon, J. Chromatogr. A, 736 (1996) 61-75.
49
50 414 [8] B. Tienpont, F. David, A. Pereira, P. Sandra, J. Chromatogr. A, 1218 (2011) 8303-8312.
51
52 415 [9] G. Montaudo, R. P. Lattimer (Eds.), Mass Spectrometry of Polymers, CRC Press, 2002.
53
54
55 416 [10] R. B. Cole (Ed.), Electrospray and MALDI Mass Spectrometry, Wiley, 2010.
56
57 417 [11] C. Barner-Kowollik, T. Gruending, J. Falkenhagen, S. Weidner (Eds), Mass Spectrometry in Polymer
58
59 418 Chemistry, Wiley-VCH, 2012.

1 419 ~~[12] T. Usami, F. Keitoku, H. Ohtani, S. Tsuge, Polymer, 33 (1992) 3024-303.~~
2
3 420 ~~[13] R. S. Lehrle, E. J. Place, Polym. Degrad. Stab., 56 (1997) 215-219; 221-226.~~
4
5 421 ~~[14] F. Bertini, G. Audisio, V. V. Zuev, Polym. Degrad. Stab, 89 (2005) 233-239.~~
6
7 422 ~~[15] H. Ito, K. Hayakawa, A. Yamamoto, A. Murase, K. Hayakawa, M. Kuno, Y. Inoue, J. Chromatogr. A,~~
8
9
10 423 ~~1132 (2006) 256-259.~~
11
12 424 ~~[16] T. Tsuge, Y. Iiguni, H. Ohtani, H. Ito, M. Kuno, Bunseki Kagaku, 59 (2010) 693-700.~~
13
14 425 ~~[17] H. Ohtani, T. Nagaya, Y. Sugimura, S. Tsuge, J. Anal. Appl. Pyrolysis, 4 (1982) 117-131.~~
15
16
17 426 [12] H. Ohtani, T. Kimura, S. Tsuge, Anal. Sci., 2 (1986) 179-182.
18
19 427 [13] T. Koshibuka, T. Ohkawa, K. Takeda, Polym. Degrad. Stab., 79 (2003) 1-11.
20
21
22 428 [14] F. Samperi, C. Puglisi, R. Alicata, G. Montaudo, Polym. Degrad. Stab., 83 (2004) 11-17.
23
24 429 [15] J. Jachowicz, M. Kryszewski, P. Kowalski, J. Polym. Sci., 22 (1978) 2891-2899.
25
26 430 [16] T. Usami, F. Keitoku, H. Ohtani, S. Tsuge, Polymer, 33 (1992) 3024-303.
27
28
29 431 ~~[23] S. J. Wetzel, C. M. Guttman, J. E. Girard, Int. J. Mass Spectrom., 238 (2004) 215-225.~~
30
31 432 [17] H. Matsubara, A. Yoshida, Y. Kondo, S. Tsuge, H. Ohtani, Macromolecules, 36 (2003) 4750-4755.
32
33 433
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65

434 Figure captions:

435

436 **Fig. 1.** Schematic diagram and operation of pyrolysis unit with movable trapping tube for pyrolyzates. (a)

437 Pyrolysis and trapping of pyrolyzates. (b) Collection of pyrolyzates by solvent flushing. (1) heating block,

438 (2) carrier gas, (3) seal gas, (4) cooling fan, (5) movable pyrolysis tube, (6) sample cup, (7) flushing solvent,

439 (8) 3-port switching valve, (9) Kalrez O-ring (10) quartz wool, (11) solvent collection vial, (12) PTFE tube,

440 (13) stainless steel filter, (14) gas outlet, and (15) PTFE 3-port valve.

441

442 ~~Fig. 2. Operation of the pyrolysis unit. (a) Pyrolysis and trapping of pyrolyzates. (b) Collection of~~

443 ~~pyrolyzates by solvent flushing.~~

444

445 ~~Fig. 3. Polymer samples.~~

446

447 ~~Fig. 4. Pyrogram of PMA observed by Py-GC-MS at 550 °C.~~

448

449 ~~Fig. 5. MALDI mass spectrum of pyrolysis products at 550 °C of PMA collected in THF. (a) overall~~

450 ~~spectrum. (b) its expended view between $m/z = 1550$ and 1600 .~~

451

452 ~~Fig. 6. Pyrogram of Nylon11 observed by Py-GC-MS at 550 °C.~~

453

454 ~~Fig. 7. MALDI mass spectrum of pyrolysis products at 550 °C of Nylon11 collected in THF. (a) overall~~

455 ~~spectrum. (b) its expended view between $m/z = 620$ and 645 .~~

456

457 ~~Fig. 8. Pyrolysis pathway of PBT.~~

458

459 **Fig. 2.** Pyrogram of PBT observed by Py-GC-MS at 500 °C.

1 460
2
3 461
4
5 462
6
7 463
8
9
10 464
11
12 465
13
14 466
15
16
17 467
18
19 468
20
21
22 469
23
24 470
25
26
27 471
28
29 472
30
31 473
32
33
34 474
35
36 475
37
38 476
39
40
41 477
42
43 478
44
45 479
46
47
48 480
49
50 481
51
52
53 482
54
55 483
56
57 484
58
59
60
61
62
63
64
65

Fig. 3. MALDI mass spectrum of pyrolysis products at 500 °C of PBT collected with THF.

Fig. 4. On-line Py-HPLC-UV chromatogram of collected pyrolysis products of PBT at 500 °C.

Fig. 5. Typical pyrolysis pathway of PPE .

Fig. 6. Pyrogram of PPE observed by Py-GC-MS at 550 °C.

Fig. 7. MALDI mass spectrum of pyrolysis products at 550 °C of PPE collected in THF.

~~**Fig. 15.** Pyrolysis pathway of PPE to form pentamer (n = 5) with maximum numbers of substitution (5-0) through complete Fries rearrangement.~~

Fig. 8. On-line HPLC chromatogram of collected pyrolysis products of PPE at 550 °C by UV detection.

~~**Fig. 17.** Total ion chromatogram of collected pyrolysis products of PPE at 550 °C obtained by HPLC-ESI-MS.~~

~~**Fig. 18.** Accumulated ESI mass spectrum of collected pyrolysis products of PPE at 550 °C. (Partial view around pentamer) –~~

Fig. 9. Extracted ion chromatograms of oligomeric pyrolysis products of PPE at 550 °C for various degrees of polymerization (n = 4-6) ~~with maximum numbers of substitution (n = 0).~~

485 **Fig. 10.** Thermally-assisted hydrolysis and methylation pathway of photocured acrylic resin sample.

1
2 486

3 487 **Fig. 11.** MALDI mass spectrum of THM products of photocured acrylic resin sample collected in THF.

4
5 488

6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65

489

Table 1

490

Representative pyrolyzates of Nylon 11.

491

peak code	chemical structure
L_n	$\left[\left(\text{CH}_2 \right)_{10} \text{NH} \right]_n$
$A_{X-(m+n)}$	$\text{CH}_3 \left(\text{CH}_2 \right)_{m-1} \text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{NH} \right]_X \left(\text{CH}_2 \right)_{n-1} \text{CH}_3$
$B_{X-(m+n)}$	$\text{H}_2\text{C}=\text{CH} \left(\text{CH}_2 \right)_{m-2} \text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{NH} \right]_X \left(\text{CH}_2 \right)_{n-1} \text{CH}_3$
$C_{X-(m+n)}$	$\text{H}_2\text{C}=\text{CH} \left(\text{CH}_2 \right)_{m-2} \text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{NH} \right]_X \left(\text{CH}_2 \right)_{n-2} \text{CH}=\text{CH}_2$
$m+n \leq 20, m \leq 10, n \leq 10$	
D_{X-n}	$\text{CH}_3 \left(\text{CH}_2 \right)_{n-1} \text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{NH} \right]_X \left(\text{CH}_2 \right)_{10} \text{CN}$
E_{X-n}	$\text{H}_2\text{C}=\text{CH} \left(\text{CH}_2 \right)_{n-1} \text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{NH} \right]_X \left(\text{CH}_2 \right)_{10} \text{CN}$
$n \leq 10$	

501

502

1 503
 2
 3 504
 4
 5 505
 6
 7 506
 8
 9 507
 10
 11 508
 12
 13 509
 14
 15 510
 16
 17 511
 18
 19 512
 20
 21 513
 22
 23
 24
 25
 26
 27
 28
 29
 30
 31
 32
 33
 34
 35
 36
 37
 38
 39
 40
 41
 42
 43
 44
 45
 46
 47
 48
 49
 50
 51
 52
 53
 54
 55
 56
 57
 58
 59
 60
 61
 62
 63
 64
 65

Table 1
Possible pyrolyzates of PBT.

peak code	chemical structure	MW
A _n		176+220n
B _n		220+220n
C _n		274+220n
D _n		166+220n
E _n		122+220n
F _n		298+220n

514

Table 2


515

Peak assignment of PPE pyrolyzates observed in dimer and trimer regions in the pyrogram.


516

chemical structure


517



MW=242




MW=256




MW=228

518


n=2



MW=242



MW=242




MW=256


519

520

521



MW=362




MW=376

522


523

524

n=3




MW=348




MW=362

525

526



MW=362



MW=376

527

528

529
1
2
3 530
4
5 531
6
7 532
8
9
10 533
11
12 534
13
14 535
15
16 536
17
18
19 537
20
21 538
22
23
24 539
25
26 540
27
28
29 541
30
31 542
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65

Table 3

Possible isomeric pyrolyzates for pentamer (5-0) of PPE observed in extracted ion chromatogram at m/z 639.

polarity	Structure of the pyrolyzates of PPE
high	
middle	
low	

1 A novel analytical pyrolysis device applicable for measurements of less volatile pyrolyzates

2

3 Hirohisa Kano^a, Takahiro Okamoto^a, Shinya Kitagawa^a, Yoshinori Iiguni^a, Hajime Ohtani^{a,*}, Hiroshi Ito^b,

4 Koichiro Iwai^{a,b}, Minoru Kuno^c

5

6 ^a *Department of Materials Science and Engineering, Graduate School of Engineering, Nagoya Institute of*

7 *Technology, Nagoya 466-8555, Japan*

8 ^b *Toyota Central R&D Labs., Inc., Nagakute, Aichi, 480-1192, Japan*

9 ^c *GL Science Inc., 237-2, Sayamagahara, Iruma, Saitama, 358-0032, Japan*

10

11 ABSTRACT

12 A novel analytical pyrolysis device suitable for the measurements of the less volatile pyrolyzates was
13 developed. This device is comprised of a pyrolysis unit based on the vertical microfurnace pyrolyzer and a
14 movable pyrolysis tube made of stainless steel. The upper side of the pyrolysis tube is first set at the heated
15 center of the pyrolysis unit. A polymer sample in a sample cup is introduced into the heated center and
16 pyrolyzed under the flow of nitrogen carrier gas. The pyrolyzates are then transferred to the bottom side of
17 the pyrolysis tube which is cooled down by air blow and the less volatile pyrolyzates are condensed. After
18 that the pyrolysis tube is moved down out of the pyrolysis unit and blocked from the hot zone by moving
19 the three port valve. Finally the trapped pyrolyzates are flushed by an appropriate solvent into a collection
20 vial. The recovered pyrolyzates in the solution can be then subjected to the measurements suitable for less
21 volatile products such as matrix-assisted laser desorption/ionization mass spectrometry. Moreover, this
22 device is also on-line connected to the injection port of a high performance liquid chromatographic system,
23 which is also appropriate for the measurements of oligomeric products. The basic performance of the
24 developed device was tested by using some polymer samples forming less volatile pyrolyzates such as
25 poly(butylene terephthalate) and poly(2,6-dimethyl-1,4-phenylene ether). In addition, the usefulness of this
26 method was demonstrated through the characterization of the crosslinking sequences in a photocured resin

1 27 combined with the thermally-assisted hydrolysis and methylation reaction.

2
3 28

4
5 29 *Keywords:* Less volatile pyrolyzate, HPLC, MALDI-MS, Pyrolysis device, Analytical pyrolysis, Polymer
6
7 30 characterization

8
9
10 31

11 32 **1. Introduction**

12
13
14 33 Recently, modern pyrolysis-gas chromatography (Py-GC) is widely utilized in the field of polymer
15 34 characterization not only for simple “fingerprint” identification and compositional analysis, but also for the
16
17 35 detailed structural analysis such as copolymer sequences, region- and stereoregularities, end groups, and
18
19 36 branching and cross-linking structures [1-3]. For such advanced objectives, it is often demanded to inspect
20
21 37 a wide range of pyrolyzates including large molar mass (oligomeric) products and/or highly polar
22
23 38 compounds. This technique, however, is inherently inapplicable for the measurements of such less volatile
24
25 39 pyrolyzates due to the use of GC separation, even if they possess useful information for the
26
27 40 characterization of the sample materials. Therefore, an analytical pyrolysis technique suitable for such
28
29 41 intractable pyrolyzates has been required.

30
31
32
33
34 42 A possible approach to overcome this issue is the combination of pyrolysis and high performance
35
36 43 liquid chromatographic (HPLC) separation. Stepwise (off-line) analyses of collected pyrolyzates by
37
38 44 solvent-trap or extraction have been reported by using HPLC (or supercritical fluid chromatography (SFC))
39
40 45 for polymeric samples such as polystyrene [4, 5] and latex containing polar vinyl acid components [6]. In
41
42 46 1996 van der Hage and Boon [7] proposed an on-line Curie-point pyrolysis-HPLC-mass spectrometry (MS)
43
44 47 system. In this method, by controlled solvent flow with appropriate valve switching, condensed pyrolyzates
45
46 48 in the pyrolysis chamber were washed out into a precolumn for trapping. The collected products were then
47
48 49 flushed into the separation column of HPLC by switching the solvent and its flowing path with the valve
49
50 50 control. The feasibility and the potential of this technique were demonstrated through the analysis of high
51
52 51 molar mass pyrolysis products of lignin polymers. Recently, a generic pyrolysis interface called “PyroVial”
53
54 52 was presented [8]. The sample in the closed vial was electrically heated up to 1000 °C on the dedicated

53 socket. Combined with a sophisticated syringe manipulation, the vial containing the pyrolyzates can be
54 hyphenated to any chromatographic methods including GC, HPLC and SFC.

55 Alternative measure for this subject would be the direct combination between pyrolysis and MS
56 without separation. In the conventional pyrolysis-MS technique, however, vaporized pyrolyzates have been
57 predominantly measured using electron ionization (EI), chemical ionization (CI) and field ionization (FI)
58 sources [9]. Less volatile pyrolyzates, therefore, cannot be analyzed and considerable fragmentation
59 sometimes makes the observed mass spectra too complicate to interpret especially for EI source.
60 Meanwhile, matrix assisted laser desorption/ionization (MALDI)-MS and electrospray ionization (ESI)-MS
61 are widely utilized in recent years for the characterization of large molar mass compounds especially for
62 oligomeric component [10, 11]. Thus the combination between pyrolysis and MALDI- and/or ESI-MS is
63 also a promising approach to obtain the information of less volatile pyrolyzates.

64 In this work, a novel pyrolysis device suitable for the analysis of the less volatile pyrolyzates was
65 developed on the basis of a microfurnace pyrolyzer which is one of the representative pyrolyzers used for
66 analytical purpose [3]. Pyrolysis reaction is achieved basically under the same conditions for conventional
67 Py-GC measurements. The device can trap the less volatile pyrolyzates in a movable pyrolysis tube and
68 flush the trapped components with an appropriate solvent. The recovered pyrolyzates in the solvent are
69 readily analyzed by MALDI- and ESI-MS. Furthermore, this device can be also on-line connected to the
70 injection port of HPLC. The basic performance and applicability of the device was tested by using some
71 polymer samples forming less volatile pyrolyzates.

72 73 **2. Experimental**

74 *2.1. Construction and operation of pyrolysis device*

75 Fig. 1 shows the schematic diagram and the operation of the developed pyrolysis device. This device
76 is comprised of a pyrolysis unit based on a vertical microfurnace pyrolyzer [3], of which temperature is
77 controlled by a PID controller (EC5530, Ohkura) and a movable pyrolysis tube (5) made of stainless steel
78 (190 mm long x 4 mm o.d. x 2 mm i.d.) set inside of an outer tube (1/4 inch o.d. x 4.35 mm i.d.). The upper

1 79 side of the pyrolysis tube is first positioned at the heated center of the pyrolysis unit (Fig. 2(a)). A polymer
2
3 80 sample weighed about 1 mg in a surface deactivated stainless steel sample cup (6) (Eco cup, Frontier Lab)
4
5 81 is introduced into the heated center at 500 or 550 °C and pyrolyzed under the flow of nitrogen carrier gas
6
7 82 (70 mL min⁻¹). As for thermally-assisted hydrolysis and methylation (THM) of a photocured resin sample,
8
9
10 83 3 mg of the resin mixed with 5 μL of a 25 wt% of tetramethylammonium hydroxide (TMAH) (Aldrich) in
11
12 84 the sample cup was introduced into the pyrolysis unit set at 400 °C. The pyrolyzates are then transferred to
13
14 85 the bottom side of the pyrolysis tube which is cooled down by cooling fans (4) and the less volatile
15
16
17 86 pyrolyzates are condensed. To increase the trapping efficiency, the connecting union at the bottom of the
18
19 87 pyrolysis tube is packed with quartz wool (10) (3 mg). During this step nitrogen seal gas (10 mL min⁻¹) (3)
20
21
22 88 flows to prevent the leakage of pyrolysis products outside of the pyrolysis (inner) tube. After that the
23
24 89 pyrolysis tube is moved down out of the pyrolysis unit sealed with a tough O-ring (9) (Kalrez S004, 3.5
25
26 90 mm x 1.5 mm, DuPont) and blocked from the hot zone by switching the three port valve (8) (Fig. 2(b)).
27
28
29 91
30

31 92 **Fig. 1.** Schematic diagram and operation of pyrolysis unit with movable trapping tube for pyrolyzates.
32
33 93

34
35 94 Finally the trapped pyrolyzates are flushed by an appropriate solvent (7) into a collection vial (5 mL
36
37 95 volume) (11). The flushing solvent (1 mL THF; HPLC grade; Tokyo Chemical Industry) is pushed out to
38
39
40 96 the pyrolysis tube with the seal gas by rotating another 3-port switching valve (8) for the seal gas. The
41
42 97 recovered pyrolyzates in the solution are readily subjected to the measurements for less volatile products
43
44
45 98 such as MALDI- and ESI-MS. Moreover, this device can be also on-line connected to the injection port of
46
47 99 HPLC. After the PTFE 3-port valve (15) is switched and the position of the collection vial is adjusted, the
48
49 100 solution of the recovered pyrolyzates is transferred to the sample loop (2 μL volume) of the injector
50
51
52 101 (Rheodyne 7725i) with the pressure of the seal gas.
53

54 102

55 103 2.2. Conditions for on-line HPLC

56

57 104 On-line HPLC measurements of the collected pyrolyzates were carried out with an LC-VP HPLC
60
61
62
63
64
65

1 105 system (Shimadzu) consisting of dual solvent delivery pumps (LC-10ADvp; Shimadzu), a degasser
2
3 106 (DGU-12A, Shimadzu), a column oven (CO-10ADvp, Shimadzu), a photodiode array UV-visible detector
4
5 107 (SPD-M10A, Shimadzu), and a system controller (SCL-10A, Shimadzu). The separation column was
6
7 108 Chemcobond 3.5-ODS-H (Chemco, 150 mm x 1.0 mm, 3.5 μm particles) connected with a guard column
8
9
10 109 of the same packings (Chemco, 35 mm x 1.0 mm). A mixed solvent of water and THF was adopted as the
11
12 110 mobile phase for a gradient elution. The gradient conditions were programmed from 30 (for PBT) or 20%
13
14 111 (for PPE) THF (10 min hold) to 100% THF (for 60 min), occasionally added with 0.1% (v/v) trifluoroacetic
15
16
17 112 acid mainly for polar products from PBT. The column oven was set at 40 $^{\circ}\text{C}$. The UV detection was carried
18
19 113 out either at 240 (for PBT) or 280 (for PPE) nm.
20
21
22 114

24 115 2.3. Conditions for MALDI-MS

26 116 In this work, either high resolution spiral time-of-flight (TOF)-MS (JMS-S3000, JEOL) or
27
28
29 117 conventional TOF-MS (AXIMA-CFRplus, Shimadzu Kratos) was used. The spiral TOF-MS was equipped
30
31 118 with Nd:YLF laser (349 nm) and performed with 0.5 ns pulse width and 250 Hz frequency, while the
32
33
34 119 conventional TOF-MS was with nitrogen laser (337 nm), operated in 3 ns pulse width and 10 Hz frequency.
35
36 120 Ions generated by the laser desorption were introduced into the flight tube with an accelerating voltage of
37
38
39 121 20 kV and separated in a positive spiral mode (JMS-S3000) and a linear mode (AXIMA-CFR plus). The
40
41 122 mass spectra were acquired by averaging 30-100 and 100 individual laser shots, with the delay time set at
42
43 123 150 and 80 ns, respectively for the spiral and the conventional TOF. Mass calibration was accomplished
44
45 124 either by using poly(ethylene glycol) (average molecular weight 400 - 2000) as external standard or a
46
47
48 125 hindered amine light stabilizer (Chimassorb 119FL, Ciba Specialty Chemicals, m/z 2285.61 for the
49
50 126 monoisotopic ion) as the internal standard. As for the matrix reagent, 2-(4-hydroxyphenylazobenzoic acid)
51
52
53 127 (HABA; Sigma-Aldrich), 4-hydroxybenzylidene malononitrile (HBMN; Alfa Aesar) and
54
55 128 trans-2-[3-(4-tert-butylphenyl)-2-methyl-2-propenylidene]malononitrile (DCTB; Tokyo Chemical Industry)
56
57
58 129 were used. Each matrix solution (10 mg mL^{-1}) was prepared by dissolving 1 mg of the matrix reagent in 0.1
59
60 130 mL of THF (HPLC grade; Tokyo Chemical Industry), then 10 μL each of the matrix and the sample
61
62
63
64
65

1 131 solutions were mixed, and an aliquot (1 μL) of the mixed solution dropped on the sample plate was
2
3 132 subjected to MALDI-MS measurement. In the case of internal calibration, 1 μL of the calibrant solution (1
4
5 133 mg mL^{-1} in THF) was also added to the sample/matrix mixed solution. The ions observed in the MALDI
6
7 134 mass spectra were predominantly sodium-cationized molecules of the products under the conditions in this
8
9
10 135 work.

11 12 136 13 14 137 *2.4. Conditions for HPLC-ESI-MS*

15 138 In order to identify the collected pyrolyzates of poly(phenylene ether) (PPE) subjected to HPLC
16
17 139 measurements, a HPLC-MS system consisting of UPLC (Acuity, Waters) and TOF-MS (SYNAPT G2,
18
19 140 Waters) was used. MS was operated in the positive ESI mode under the conditions of source temperature
20
21 141 $120\text{ }^\circ\text{C}$, desolvation N_2 gas temperature $350\text{ }^\circ\text{C}$, desolvation gas flow 800 L h^{-1} , capillary voltage of 3.0 kV ,
22
23 142 and cone voltage 40 V . The separation column and the conditions for measurements were basically the
24
25 143 same as those for the on-line HPLC-UV detection (section 2.2).
26
27
28
29
30

31 144 32 33 145 *2.5. Conditions for Py-GC-MS*

34 146 Py-GC-MS measurements of the polymer samples were also performed for comparison. A vertical
35
36 147 microfurnace pyrolyzer (PY-2020iD, Frontier Lab) was attached to the injection port of a GC-MS system
37
38 148 (GCMS-QP2010 Plus, Shimadzu). About 0.1 mg of a polymer sample weighed in a deactivated stainless
39
40 149 steel sample cup was pyrolyzed at 500 or $550\text{ }^\circ\text{C}$ under the flow of helium carrier gas. The flow rate of 50
41
42 150 mL min^{-1} of carrier gas at pyrolyzer was reduced to 1.0 mL min^{-1} at a metal capillary separation column
43
44 151 (Ultra-Alloy 5 MS, $30\text{ m} \times 0.25\text{ mm i.d.}$, coated with $0.25\text{ }\mu\text{m}$ of 5% diphenyl-95% dimethyl-polysiloxane,
45
46 152 Frontier Lab) by means of a splitter. The temperature both of PY/GC interface and the GC injection port
47
48 153 was set at $320\text{ }^\circ\text{C}$. The column oven temperature was first held at $40\text{ }^\circ\text{C}$ for 2 min , then programmed to
49
50 154 $320\text{ }^\circ\text{C}$ at a rate of $10\text{ }^\circ\text{C min}^{-1}$ and finally held at $320\text{ }^\circ\text{C}$ for 30 min . The separated components were
51
52 155 introduced into the ionization source of a quadrupole MS through a GC/MS interface held at $320\text{ }^\circ\text{C}$. The
53
54 156 ionization was achieved with electron ionization at 70 eV . The mass scan range was $m/z\ 29 - 600$.
55
56
57
58
59
60
61
62
63
64
65

157

2.6. Polymer samples

Poly(methyl acrylate) (PMA) sample was radically polymerized in a standard procedure. Industrially available pelletized samples were used for Nylon 11 (polyundecanoamide) and poly(butylene terephthalate) (PBT). Poly(2,6-dimethyl-1,4-phenylene ether) (PPE) sample obtained using Mn-based catalyst was the same as that used in previous work [12]. The results for PMA and Nylon 11 are discussed as supplementary information.

The photocured acrylic resin sample was prepared with neopentyl glycol diacrylate through photopolymerization in the presence of 7 wt% 2,2'-dimethoxy-2-phenylacetophenone as a photoinitiator using the UV-cured equipment consisting of a conveyor and a medium-pressure mercury lamp.

3. Results and discussion

3.1 Basic performance of the developed pyrolysis device

3.1.1. Poly(butylene terephthalate) (PBT)

The basic performance of the developed pyrolysis device including on-line HPLC measurements of the pyrolyzates obtained by this device was evaluated using PBT, which is widely utilized as one of typical engineering plastics and contains aromatic moieties suitable for UV detection. In general, polyesters like PBT thermally decompose into oligomeric products primarily through *cis*-elimination at the ester linkages followed by partial decarboxylation. It has been reported for PBT, therefore, several types of products listed in Table 1 should be formed by pyrolysis reactions [12-14].

Table 1 Possible pyrolyzates of PBT.

Fig. 2 shows a pyrogram of PBT obtained by Py-GC-MS at 500 °C. The main pyrolyzates observed are 1,3-butadiene and A_n, B_n and C_n with degree of polymerization (n) of 0 or 1 [3, 12]. The yields of E_n and F_n products having phenyl terminal(s) were fairly small, probably because the contribution of

183 decarboxylation to form phenyl terminals would be minor under the conditions of relatively low pyrolysis
184 temperature at around 500 °C. Meanwhile highly polar decarboxylic acid products (D_n) only shows a broad
185 hump around 20 min for terephthalic acid (D_0). At least the larger pyrolyzates than C_1 , or the products with
186 molar mass > 500, were not able to be observed by conventional Py-GC-MS at all.

Fig. 2. Pyrogram of PBT observed by Py-GC-MS at 500 °C.

190 Next PBT was pyrolyzed at 500 °C in the pyrolysis device, and the pyrolyzates were collected in THF.
191 The obtained THF solution was first measured directly by conventional MALDI-MS with HBMN as a
192 matrix reagent. Fig. 3 shows the MALDI mass spectrum of the recovered pyrolyzates of PBT. The
193 pyrolyzates observed in the pyrogram by Py-GC-MS (Fig. 2) such as B_0 , A_1 - C_1 and E_1 were also observed,
194 although smaller products (A_0 , E_0 , etc) might evaporate during MALDI process and not be detected in the
195 mass spectrum. In addition, the pyrolyzates that could be hardly detected by Py-GC-MS, *e. g.* D_1 , F_1 and
196 the products with greater degrees of polymerization ($n = 2, 3$), were also clearly confirmed in the MALDI
197 mass spectrum (Fig. 3). These observations prove that such larger pyrolyzates from PBT can be sufficiently
198 recovered by the pyrolysis device.

Fig. 3. MALDI mass spectrum of pyrolysis products at 500 °C of PBT collected with THF.

202 The solution of the collected pyrolyzates was then subjected to on-line HPLC measurements. Fig. 4
203 shows the UV-detected chromatogram of the recovered pyrolyzates of PBT. In this chromatogram, several
204 peaks were clearly observed. In order to identify the peak components, the eluents around the main peaks
205 (26 - 29, 34 - 37, 42 - 44, and 46 - 48 min) were fractionated (5 times for accumulation) and subjected to
206 MALDI-MS measurements. As the results, the pyrolysis products designated in Fig. 4 were at least
207 assigned. Among these, the highly polar and/or large molar mass components, which cannot be detected by
208 Py-GC-MS, were included (A_2 , B_2 , C_2 , D_2 , E_2 , etc). This result demonstrates that on-line HPLC

209 measurements of less volatile pyrolyzates are able to be performed by using the analytical pyrolysis system
developed.

Fig. 4. On-line Py-HPLC-UV chromatogram of collected pyrolysis products of PBT at 500 °C.

3.1.2. Poly(2,6-dimethyl-1,4-phenylene ether) (PPE)

The pyrolysis device was then applied to the characterization of PPE. It has been reported that PPE thermally decompose into various oligomers in a unique manner with the major contribution of Fries rearrangement reaction to form C-C bondings along with the simple cleavages of the main chains as shown in Fig. 5 [15, 16]. Fig. 6 shows a typical pyrogram of PPE obtained by Py-GC-MS at 550 °C. Although a wide range of oligomeric products should be formed, observed pyrolyzates by Py-GC-MS were limited to mainly monomers, dimers and trimers together with minor amounts of tetramers. Table 2 summarizes the peak assignments of major pyrolyzates in dimeric and trimeric regions. Reflecting the considerable contribution of Fries arrangement, various types of products containing some isomers with methylene and/or ether linkage(s) can be formed [3, 15, 16].

Fig. 5. Typical pyrolysis pathway of PPE .

Fig. 6. Pyrogram of PPE observed by Py-GC-MS at 550 °C.

Table 2 Peak assignment of PPE pyrolyzates observed in dimer and trimer regions in the pyrogram.

The PPE sample was also pyrolyzed in the new pyrolysis devices at 550 °C. The recovered solution of the pyrolyzates in THF was first analyzed by MALDI-spiral-TOF-MS using DCTB as the matrix reagent. The observed mass spectrum is shown in Fig. 7. In addition to the pyrolyzates up to tetramers (n = 4) observed by Py-GC-MS, higher molar mass oligomers are clearly detected at least up to heptamers (n = 7). This fact shows the less volatile pyrolyzates of PPE can be successfully recovered by the pyrolysis device.

Fig. 7. MALDI mass spectrum of pyrolysis products at 550 °C of PPE collected in THF.

236

237 The recovered solution of the PPE pyrolyzates was then subjected to on-line HPLC measurements. Fig.
238 8 shows the chromatogram of the collected PPE pyrolyzates observed in UV detection at 280 nm. The fact
239 that a number of separated peaks were clearly observed suggest that the analysis of the less volatile
240 pyrolyzates for PPE can be achieved by on-line HPLC measurement. In the case of the PPE pyrolyzates,
241 however, the correct identification of each eluent should be difficult by the fractionation-MALDI-MS
242 technique because of possible decomposition and secondary reaction during laser irradiation as mentioned
243 in supplementary information. Therefore, HPLC-ESI-MS measurements were carried out for the
244 appropriate identification of the PPE pyrolyzates, in which an aliquot of the collected PPE pyrolyzate
245 solution was sampled and introduced separately (off-line) into HPLC-ESI-MS.

246

Fig. 8. On-line HPLC chromatograms of collected pyrolysis products of PPE at 550 °C by UV detection.

248

249 The elution behaviors of individual pyrolyzates were examined in detail based on the corresponding
250 selected ion chromatograms (SIMs) obtained by HPLC-ESI-MS. Fig. 9 shows the SIMs of some PPE
251 pyrolyzates for the largest components of individual oligomers (tetramers to hexamers). Not only the
252 elution times of the products shifted depending on their molecular size, but also the number of the observed
253 peaks for the individual oligomers increased with increase in m/z . The latter multiple peaks should be
254 attributed to the possible formation of the isomeric pyrolyzates. As for pentamers (m/z 639), at least five
255 separated peaks are observed. If the pyrolysis of PPE proceeds with the complete Fries rearrangement, the
256 formed pentamer must be that containing 5 OH groups (Fig. S5). However, considerable amounts of
257 dimeric and trimeric pyrolyzates containing ether linkage(s) were also confirmed in the pyrogram of PPE
258 obtained by Py-GC-MS as shown in Fig. 6 and Table 2. This fact suggests that various isomeric products
259 containing ether linkage(s) should be also formed in the pentamer region. Even for the pentameric products
260 with the same m/z 639, various isomeric forms can be considered as shown in Table 3. (Further isomeric
261 products should be produced depending on the difference in the distributions of OH groups.) Since these

262 pentamers contain various numbers of OH groups (1 - 5), their polarity should be different from each other
263 which in turn causes the peak separation in the chromatogram. This result demonstrated that even the
264 isomeric products in less volatile pyrolyzates can be separated by HPLC measurements which could not be
265 differentiated by MS analyses alone.

Fig. 9. Extracted ion chromatograms of oligomeric pyrolysis products of PPE at 550 °C.

Table 3 Possible isomeric pyrolyzates for pentamer of PPE observed in extracted ion chromatogram.

3.2 Application to elucidate crosslinking sequences in photocured resin using THM

Photocurable acrylic esters are utilized for many industrial materials. Characterization of the network structures in the cured resins has been requested to correlate it with various properties of the resins. It has been demonstrated that THM-GC often brings us valuable clues on the chemical structure in the photocured acrylic resins. In this technique, the ester linkages in the cured resin samples selectively decomposed into methyl derivatives which provide various information such as conversion, copolymer composition and average molar mass of prepolymers [2]. As for the cross-linking sequences, it was reported that the chain length distribution of the network junctions composed of up to six acryloyl units in the photocured acrylic resin sample could be estimated from the peak intensities of the characteristic methyl acrylate (MA) oligomers (up to hexamer) reflecting cross-linking sequence structure [17]. However, the network junctions containing longer sequence comprised of more than seven acryloyl units were not elucidated because the corresponding higher MA oligomers were not observed in the pyrogram due to their low volatility.

In this case, the developed pyrolysis device could provide the information of the longer cross-linking sequences in the cured acrylic resins through the THM reaction in the presence of TMAH. Fig. 10 presents possible THM reactions in the photocured resin sample formed from neopentylglycol diacrylate, in which two types of MA oligomers might be produced. The cured resin sample mixed with TMAH was introduced into the pyrolysis device at 400 °C. The collected products in THF were then analyzed by conventional MALDI-MS with HABA as a matrix reagent and NaI as a cationization agent. Fig. 11 shows the MALDI

1 288 mass spectrum of the THM products of the cured acrylic resin sample. Two series of the MA oligomers
2
3 289 were observed at least up to 18-mers. This result demonstrates that the developed pyrolysis device can
4
5 290 provide us the information of the chain length distribution of the network junctions in cured resins even for
6
7 291 the longer sequences which cannot be evaluated by THM-GC measurements.
8
9

10 292
11
12 293 **Fig. 10.** Thermally-assisted hydrolysis and methylation pathway of photo cured acrylic resin sample.
13

14 294 **Fig. 11.** MALDI mass spectrum of THM products of photocured acrylic resin sample collected in THF.
15
16

17 295 18 19 296 **4. Conclusion**

20
21
22 297 A novel pyrolysis device suitable for analysis of the less volatile pyrolyzates was developed using a
23
24 298 microfurnace type pyrolysis unit combined with a movable pyrolysis tube (made of stainless steel). The
25
26 299 less volatile pyrolyzates of polymer samples were trapped on the inner wall of the pyrolysis tube by air
27
28 300 cooling, and then solvent-flushed into a collection vial. The recovered solution containing less volatile
29
30 301 pyrolyzates can be directly analyzed by MALDI-MS. The on-line HPLC analyses of the solutions were also
31
32 302 possible. In the present state, however, only the HPLC-UV detection system has been on-line connected
33
34 303 with the pyrolysis device. The more powerful combination could be that between pyrolysis and HPLC-MS.
35
36 304 The development of on-line connection of HPLC-ESI-MS with the pyrolysis device is currently in progress
37
38 305 along with the optimization of the operating conditions.
39
40
41 306
42
43
44

45 307 **Acknowledgments**

46 308 This work was supported in part by the research foundation of Frontier Laboratory Inc., and JPSJ
47
48 309 KAKENHI Grant Number 233350032.
49
50
51
52
53 310

54 311 **References**

- 55
56
57 312 [1] S. C. Moldveanu, *Analytical Pyrolysis of Synthetic Organic Polymers*, Elsevier, 2005.
58
59 313 [2] T. P. Wampler (Ed.), *Applied Pyrolysis Handbook*, second ed., CRC Press, 2007.
60

- 1 314 [3] S. Tsuge, H. Ohtani, C. Watanabe, Pyrolysis-GC/MS Data Book of Synthetic Polymers -Pyrograms,
2
3 315 Thermograms and MS of Pyrolyzates-, Elsevier, 2011.
4
5 316 [4] S.-T. Lai, D.C. Locke, J. Chromatogr. A, 255 (1983) 511-527.
6
7 317 [5] T. Nonobe, H. Ohtani, T. Usami, T. Mori, H. Fukumori, Y. Hirata, S. Tsuge, J. Anal. Appl. Pyrolysis, 33
8
9 318 (1995) 121-138.
10
11
12 319 [6] F. C.-Y. Wang, B. Gerhart, C. G. Smith, Anal. Chem., 67 (1995) 3681-3686.
13
14
15 320 [7] E.R.E. van der Hage, J.J. Boon, J. Chromatogr. A, 736 (1996) 61-75.
16
17 321 [8] B. Tienpont, F. David, A. Pereira, P. Sandra, J. Chromatogr. A, 1218 (2011) 8303-8312.
18
19 322 [9] G. Montaudo, R. P. Lattimer (Eds.), Mass Spectrometry of Polymers, CRC Press, 2002.
20
21
22 323 [10] R. B. Cole (Ed.), Electrospray and MALDI Mass Spectrometry, Wiley, 2010.
23
24 324 [11] C. Barner-Kowollik, T. Gruending, J. Falkenhagen, S. Weidner (Eds), Mass Spectrometry in Polymer
25
26 325 Chemistry, Wiley-VCH, 2012.
27
28
29 326 [12] H. Ohtani, T. Kimura, S. Tsuge, Anal. Sci., 2 (1986) 179-182.
30
31 327 [13] T. Koshibuka, T. Ohkawa, K. Takeda, Polym. Degrad. Stab., 79 (2003) 1-11.
32
33
34 328 [14] F. Samperi, C. Puglisi, R. Alicate, G. Montaudo, Polym. Degrad. Stab., 83 (2004) 11-17.
35
36 329 [15] J. Jachowicz, M. Kryszewski, P. Kowalski, J. Polym. Sci., 22 (1978) 2891-2899.
37
38 330 [16] T. Usami, F. Keitoku, H. Ohtani, S. Tsuge, Polymer, 33 (1992) 3024-303.
39
40
41 331 [17] H. Matsubara, A. Yoshida, Y. Kondo, S. Tsuge, H. Ohtani, Macromolecules, 36 (2003) 4750-4755.
42
43 332
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65

1 333 Figure captions:
2
3 334
4
5 335 **Fig. 1.** Schematic diagram and operation of pyrolysis unit with movable trapping tube for pyrolyzates. (a)
6
7 336 Pyrolysis and trapping of pyrolyzates. (b) Collection of pyrolyzates by solvent flushing. (1) heating block,
8
9
10 337 (2) carrier gas, (3) seal gas, (4) cooling fan, (5) movable pyrolysis tube, (6) sample cup, (7) flushing solvent,
11
12 338 (8) 3-port switching valve, (9) Kalrez O-ring (10) quartz wool, (11) solvent collection vial, (12) PTFE tube,
13
14
15 339 (13) stainless steel filter, (14) gas outlet, and (15) PTFE 3-port valve.
16
17 340
18
19 341 **Fig. 2.** Pyrogram of PBT observed by Py-GC-MS at 500 °C.
20
21
22 342
23
24 343 **Fig. 3.** MALDI mass spectrum of pyrolysis products at 500 °C of PBT collected with THF.
25
26
27 344
28
29 345 **Fig. 4.** On-line Py-HPLC-UV chromatogram of collected pyrolysis products of PBT at 500 °C.
30
31 346
32
33
34 347 **Fig. 5.** Typical pyrolysis pathway of PPE .
35
36 348
37
38 349 **Fig. 6.** Pyrogram of PPE observed by Py-GC-MS at 550 °C.
39
40
41 350
42
43 351 **Fig. 7.** MALDI mass spectrum of pyrolysis products at 550 °C of PPE collected in THF.
44
45
46 352
47
48 353 **Fig. 8.** On-line HPLC chromatogram of collected pyrolysis products of PPE at 550 °C by UV detection.
49
50 354
51
52
53 355 **Fig. 9.** Extracted ion chromatograms of oligomeric pyrolysis products of PPE at 550 °C for various degrees
54
55 356 of polymerization (n = 4-6).
56
57 357
58
59
60
61
62
63
64
65

358 **Fig. 10.** Thermally-assisted hydrolysis and methylation pathway of photocured acrylic resin sample.

1
2 359

3 360 **Fig. 11.** MALDI mass spectrum of THM products of photocured acrylic resin sample collected in THF.

4
5 361

6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65

1 362
 2
 3 363
 4
 5 364
 6
 7 365
 8
 9 366
 10
 11 367
 12
 13 368
 14
 15 369
 16
 17 370
 18
 19 371
 20
 21 372
 22
 23
 24
 25
 26
 27
 28
 29
 30
 31
 32
 33
 34
 35
 36
 37
 38
 39
 40
 41
 42
 43
 44
 45
 46
 47
 48
 49
 50
 51
 52
 53
 54
 55
 56
 57
 58
 59
 60
 61
 62
 63
 64
 65

Table 1
Possible pyrolyzates of PBT.

peak code	chemical structure	MW
A _n		176+220n
B _n		220+220n
C _n		274+220n
D _n		166+220n
E _n		122+220n
F _n		298+220n

373

Table 2


374

Peak assignment of PPE pyrolyzates observed in dimer and trimer regions in the pyrogram.


375

chemical structure


376



MW=242




MW=256




MW=228

377


n=2



MW=242




MW=242




MW=256

379

380




MW=362




MW=376

381

382




MW=348




MW=362

383

n=3



MW=362



MW=376

384

385

386

387

388
1
2
389
3
4
5 390
6
7 391
8
9
10 392
11
12 393
13
14 394
15
16 395
17
18 396
19
20 397
21
22 398
23
24 399
25
26 400
27
28 401
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60
61
62
63
64
65

Table 3

Possible isomeric pyrolyzates for pentamer of PPE observed in extracted ion chromatogram at m/z 639.

polarity	Structure of the pyrolyzates of PPE
high	
middle	
low	

Fig. 1




Fig. 1. Schematic diagram and operation of pyrolysis unit with movable trapping tube for pyrolyzates

Fig. 2




Fig. 2. Pyrogram of PBT observed by Py-GC-MS at 500 °C.

Fig. 3




Fig. 3. MALDI mass spectrum of pyrolysis products at 500 °C of PBT collected with THF.

Fig. 4




Fig. 4. On-line Py-HPLC-UV chromatogram of collected pyrolysis products of PBT at 500 °C.

Fig. 5




Fig. 5. Typical pyrolysis pathway of PPE .

Fig. 6




Fig. 6. Pyrogram of PPE observed by Py-GC-MS at 550 °C

Fig. 7




Fig. 7. MALDI mass spectrum of pyrolysis products at 550 °C of PPE collected in THF.

Fig. 8




Fig. 8. On-line HPLC chromatogram of collected pyrolysis products of PPE at 550 °C by UV detection.

Fig. 9




Fig. 9. Extracted ion chromatograms of oligomeric pyrolysis products of PPE at 550 °C.

Fig. 10




Fig. 10. Thermally-assisted hydrolysis and methylation pathway of photocured acrylic resin sample.

Fig. 11




Fig. 11. MALDI mass spectrum of THM products of photocured acrylic resin sample collected in THF.

S1. Results for poly(methyl acrylate) (PMA)

PMA generally decomposes into a series of MA oligomers through the random chain scission accompanied by radical transfer reactions at elevated temperature [S1, S2]. Fig. S1 shows a typical pyrogram of PMA observed by conventional Py-GC-MS at 550 °C. In this case MA monomer and oligomers up to heptamer ($n = 7$) can be observed [S3].




Fig. S1. Pyrogram of PMA observed by Py-GC-MS at 550 °C.

PMA sample was then subjected to the developed pyrolysis device. An aliquot of the collected pyrolyzates (at 550 °C) dissolved in THF was measured by MALDI-spiral-TOF-MS using HABA as a matrix reagent. The observed mass spectrum is shown in Fig. S2. A series of MA oligomers are observed at least up to 26-mer, and each n -mer contains three types of products as shown in Fig. S2(b). These results demonstrate that high molar mass pyrolyzates unobservable by Py-GC are successfully recovered by the pyrolysis device.

S2. Results for Nylon 11

Aliphatic polyamides (Nylon) are known to be pyrolyzed mainly into various sizes of cyclic amides (lactams) and linear products having nitrile, alkyl and/or olefinic terminals formed through *cis*-elimination reactions along with random chain scissions. As for Nylon 11, therefore, the pyrolysis products listed in Table S1 should be yielded. Fig. S3 shows a typical pyrogram of Nylon 11 obtained by Py-GC-MS at 550 °C. Although various pyrolyzates listed in Table S1 are observed [S3, S4], the assigned products are mostly those having only one or no amide group such as L_1 , L_2 and D_{0-10} . Even the largest pyrolyzate, A_{1-20} , of which molecular weight is *ca.* 508, contains merely two amide groups.




Fig. S2. MALDI mass spectrum of pyrolysis products at 550 °C of PMA collected in THF.

(a) overall spectrum. (b) its expanded view between $m/z = 1550$ and 1600.

Nylon 11 was then pyrolyzed at 550 °C in the pyrolysis device and the pyrolyzates were collected with THF. The recovered solution was sampled to be measured directly by MALDI-MS. Fig. S4 shows the observed mass spectra for the pyrolyzates of Nylon 11 obtained by MALDI-spiral-TOF-MS with HABA matrix. In the oligomer range larger than m/z 500, which cannot be detected by Py-GC-MS, a number of products are detected. Among these, trimeric and tetrameric lactams (L_3 and L_4) are distinctly observed. As shown in the expanded spectrum in Fig. S4(b), various linear products containing three amide groups were identified in the m/z range between 620 and 645. (The observed peak splitting should be caused by the overlapping of possible isotopic ions for individual pyrolyzates.) Moreover, although with small intensities, ions for the pyrolyzates were clearly observed even in the region around 1400, which corresponding to the octameric products with seven amide groups. These observations indicated that not only large molar mass but also polar pyrolyzates originated from polar polymers such as Nylons should be able to be analyzed by using the developed pyrolysis device.

Table S1 Representative pyrolyzates of Nylon 11.

peak code	chemical structure
L_n	$\left[\left(\text{CH}_2 \right)_{10} \text{C}(=\text{O})\text{NH} \right]_n$
$A_{X-(m+n)}$	$\text{CH}_3 \left(\text{CH}_2 \right)_{m-1} \text{C}(=\text{O})\text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{C}(=\text{O})\text{NH} \right]_X \left(\text{CH}_2 \right)_{n-1} \text{CH}_3$
$B_{X-(m+n)}$	$\text{H}_2\text{C}=\text{CH} \left(\text{CH}_2 \right)_{m-2} \text{C}(=\text{O})\text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{C}(=\text{O})\text{NH} \right]_X \left(\text{CH}_2 \right)_{n-1} \text{CH}_3$
$C_{X-(m+n)}$	$\text{H}_2\text{C}=\text{CH} \left(\text{CH}_2 \right)_{m-2} \text{C}(=\text{O})\text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{C}(=\text{O})\text{NH} \right]_X \left(\text{CH}_2 \right)_{n-2} \text{CH}=\text{CH}_2$ $m+n \leq 20, m \leq 10, n \leq 10$
D_{X-n}	$\text{CH}_3 \left(\text{CH}_2 \right)_{n-1} \text{C}(=\text{O})\text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{C}(=\text{O})\text{NH} \right]_X \left(\text{CH}_2 \right)_{10} \text{CN}$
E_{X-n}	$\text{H}_2\text{C}=\text{CH} \left(\text{CH}_2 \right)_{n-1} \text{C}(=\text{O})\text{NH} \left[\left(\text{CH}_2 \right)_{10} \text{C}(=\text{O})\text{NH} \right]_X \left(\text{CH}_2 \right)_{10} \text{CN}$ $n \leq 10$




Fig. S3. Pyrogram of Nylon 11 observed by Py-GC-MS at 550 °C.




Fig. S4. MALDI mass spectrum of pyrolysis products at 550 °C of Nylon11 collected in THF.
 (a) overall spectrum. (b) its expanded view between $m/z = 620$ and 645 .

S3. Secondary reactions in MALDI process of PPE pyrolyzates

A number of clusters were observed in the MALDI mass spectrum of the PPE pyrolyzates for oligomers of each degree of polymerization (n) (Fig. 7). Here the peaks in the pentamer region ($n = 5$) were examined in detail. A representative formation pathway of a pentameric ($n = 5$) product is shown in Fig. S5, provided that the Fries rearrangement would completely proceed during pyrolysis pathway. The resultant pyrolyzate shown in Fig. S5 should be possibly the largest pentamer, and is abbreviated to 5-0. Depending on the cleavage position, the smaller pentamers by one or two CH_2 unit(s) could be also produced, which are abbreviated 5-1 and 5-2, respectively. These three types of products, 5-0 - 5-2, would be thus the most probable pentameric pyrolyzates. In fact, as for the dimeric and trimeric pyrolyzates of PPE observed by Py-GC-MS, the corresponding products were confirmed to be predominant as shown in Fig. S6 and Table 3. In the observed MALDI mass spectrum of the PPE pyrolyzates, however, 5-0 - 5-2 were not the main products in pentamer region as shown in the expanded view in Fig. S6. Smaller components corresponding to 5-3 - 5-6 appeared rather abundantly. Moreover, the larger m/z ions by 2 for 5-0 - 5-6, which might be attributed to a hydrogenated ($+\text{H}_2$) or an oxidized ($+\text{O} - \text{CH}_2$) product of 5-0 - 5-6, were detected in higher intensities than those for 5-0 - 5-6. These observations suggest that the undesirable degradations and/or

secondary reactions might fairly occur for PPE pyrolyzates during MALDI process [S5].




Fig. S5. Pyrolysis pathway of PPE to form pentamer ($n = 5$) with maximum numbers of substitution (5-0) through complete Fries rearrangement.




Fig. S6. Expanded MALDI mass spectrum for pentamer region of pyrolysis products at 550 °C of PPE collected in THF

Fig. S7 shows the observed total ion chromatogram (TIC) of the collected pyrolyzates of PPE observed by HPLC-ESI-MS. Since the separation conditions have not been optimized, the observed peak separation is far from sufficient. The eluted components, however, can be identified on the basis of the accumulated ESI mass spectrum for the whole range of the chromatogram. Fig. S8 shows the expanded region for the pentamers in the accumulated ESI mass spectrum in the range between 10 and 50 min. It is interesting to note that the observed main pentamers (5-0 - 5-2) are those expected to form by pyrolysis as shown in Fig. S5. This result indicates the PPE pyrolyzates can be observed in the original form by ESI-MS, unlike MALDI-MS.




Fig. S7. Total ion chromatograms of collected pyrolysis products of PPE obtained by HPLC-ESI-MS.




Fig. S8. Accumulated ESI mass spectrum of collected pyrolysis products of PPE at 550 °C.

References

- [S1] R. S. Lehrle, E. J. Place, *Polym. Degrad. Stab.*, 56 (1997) 215-219; 221-226.
- [S2] F. Bertini, G. Audisio, V. V. Zuev, *Polym. Degrad. Stab.*, 89 (2005) 233-239.
- [S3] S. Tsuge, H. Ohtani, C. Watanabe, *Pyrolysis-GC/MS Data Book of Synthetic Polymers -Pyrograms, Thermograms and MS of Pyrolyzates-*, Elsevier, 2011.
- [S4] H. Ohtani, T. Nagaya, Y. Sugimura, S. Tsuge, *J. Anal. Appl. Pyrolysis*, 4 (1982) 117-131.
- [S5] S. J. Wetzel, C. M. Guttman, J. E. Girard, *Int. J. Mass Spectrom.*, 238 (2004) 215-225.